

Bluetooth Fingertip Pulse Oximeter JPD-500G (LED)

Features

JUMPER



SpO2
Measuring



Pulse Rate
Measuring



Multiple
Display Modes



Adjustable
Brightness



Alert
Function



Automatic
Power-Off

- Compact, elegant design
- Lightweight
- Easy to read with its 0.95" LED screen
- Built-in lithium battery and USB charging
- Includes Bluetooth 4.0 function
- Real time monitoring of vital signs with mobile app
- Automatic turn off feature to save on power consumption
- Uses Photoelectric Oxyhemoglobin Inspection technology

Specifications

JUMPER

Display	LED
Parameters	PR / SPO ₂ / PI
Accuracy	SPO ₂ - 70-99% ±2% Pulse Rate - 25-250 BPM ±2 BPM
Dimensions	58 x 32 x 32.9 mm
Weight	50.4g
Battery	AAA x2
Application Range	Adult, Adolescent
Indication	Low Battery Reminder, Value Limit Reminder



JUMPER Health App

Measure



FHR



Weight



Oximeter



Thermometer




Smart Wireless
Thermometer



TENS Therapy



Blood pressure



Multi-Device Connectivity

Trend Graph

- Fast and Accurate data acquisition that reduces human error
- Easy Data Access
- Track Vital Signs without hassle

