






FINGERTIP PULSE OXIMETER

JPD-500G(OLED)

It adopts Photoelectric Oxyhemoglobin Inspection Technology in accordance with Capacity Pulse Scanning Recording Technology, which can be used to measure human oxygen saturation & pulse rate through finger.

The device is suitable for daily care in family, hospitals, oxygen bars, community healthcare and physical care in sports, etc. It can be used before or after exercise, it is recommended to measure in a calm state.




- 
 Multiple Display Modes
- 
 Pulse Rate Measuring
- 
 SpO₂ Measuring
- 
 Brightness Adjustable
- 
 Alert Function



SUITABLE FOR ALL KINDS OF PEOPLE

Designed for sports enthusiasts, aviators, bikers, mountain climbers or people interested in measuring their SpO₂ and pulse rate.




ADJUSTABLE BRIGHTNESS SCREEN



MULTIPLE DISPLAY MODES



PLETHYSMOGRAPH & PI



WIDE-RANGE FINGER SIZES



SPECIFICATIONS

DISPLAY OLED	PARAMETERS PR / SpO₂ / PI	APPLICATION RANGE ADULT / ADOLESCENT	
DIMENSIONS 58.2 x 32.2 x 32 mm	DISPLAY DIRECTION 4 DIRECTIONS	POWER SUPPLY AAA x 2	WEIGHT 51.5g
ACCURACY SpO₂	70 - 100% ± 2 %	INDICATION REAL-TIME BATTERY STATUS	
PULSE RATE	25-250 BPM ± 2 BPM	VALUE LIMITED REMIND	

PACKAGING



- ACCESSORIES
- STORAGE BAG X 1**
 - AAA BATTERIES X 2**
 - HANGING ROPE X 1**
 - USER MANUAL X 1**
 - QUICK GUIDE X 1**

CTN SIZE(CM)	QTY/CTN	GW(KG)	FCL/20
46 X 28 X 29	100	10	680