

JUMPER


JPD-500G (OLED) Fingertip Pulse Oximeter

Features

JUMPER

Specifications

JUMPER




| | |
|-------------------|---|
| Display | OLED |
| Parameters | PR / SPO ₂ / PI |
| Accuracy | SPO ₂ - 70-99% ±2% Pulse Rate - 25-250 BPM ±2 BPM |
| Dimensions | 58.2 x 32.2 x 32 mm |
| Weight | 51.8g |
| Battery | AAA x2 |
| Application Range | Adult, Adolescent |
| Indication | Real-Time Battery Status, Value Limit Reminder |
| Display Direction | 4 Directions |


- Compact, elegant design
- Lightweight
- Easy to read with its 0.95" OLED screen
- Built-in lithium battery and USB charging
- Real time monitoring of vital signs
- Automatic turn off feature to save on power consumption
- Uses Photoelectric Oxyhemoglobin Inspection technology

Structure

JUMPER




ADJUSTABLE BRIGHTNESS SCREEN



MULTIPLE DISPLAY MODES



PLETHYSMOGRAPH & PI



WIDE-RANGE FINGER SIZES

