



## Bluetooth Infrared Thermometer JPD-FR400

### Features

JUMPER



Professional Accuracy



1 Second Reading



Hygienic Non-Touch



Smart Fever Alert

- Reading in 1 second
- Non-contact, BPA & Latex free Safe & Hygienic
- Store and Display 20 Readings in Fahrenheit or Celsius
- Smart Dual mode with Fever Alert
- Bluetooth with Smart Management App
- Uses Infrared Radiation Technology




### Specifications

JUMPER

Display	HD LCD Display
Measuring Range & Accuracy	Forehead Mode: 32.0-42.2°C ± 0.3°C 89.6-107.9°F ± 0.6°F Ear Mode: 0.0-100.0°C ± 0.2°C 32.0-212.0°F ± 0.4°F
Memory Storage	20 sets of Measurement data
Temperature units	°C/°F adjustable
Automatic shutdown function	10s ± 1s
Battery	AAAx2



### JUMPER Health App



### Multi-Device Connectivity

### Trend Graph

- Fast and Accurate data acquisition that reduces human error
- Easy Data Access
- Track Vital Signs without hassle

