







# INFRARED THERMOMETER

- Dual-mode thermometer - Measure temperature via forehead and ear.
- Reliable and accurate – Clinically tested and proven to be reliable and accurate.
- Professional medical equipment – It is FDA&CE cleared. The most recommended thermometer by doctors and pediatricians.
- Readings retrieval – Can store and recall up to 20 readings.
- Display reading in Fahrenheit or Celsius.



-  Bluetooth Connectivity
-  Professional Accuracy
-  1 Second Reading
-  Hygienic Non Touch
-  Smart Fever Alarm
-  Three Colors Backlight

## ◀◀ JPD-FR400

Intended Use: Forehead mode for all ages  
Eardrum mode for children above 1 year old  
Display: HD LCD Display  
Accuracy:  
Forehead: 32.0°C-42.2°C/ 89.6°F – 107.9 °F(±0.3°C /±0.6°F)  
Ear: 0.0°C-100.0°C/32.0°F – 212.0 °F(±0.2°C /±0.4°F)  
Battery: AAAx2  
Fever Warning: 6 short double beeps  
Memory Storage: 20 sets of measurement data  
Temperature Units: °C/°F adjustable  
Automatic Shutdown Function: 10s±1s  
Display Resolution : 0.1°C/0.1°F