



Bluetooth Blood Pressure Monitor JPD-HA120

Features

JUMPER



Display Date and Time



Wide cuff range



Automatically off



Irregular Heart Beat Indicator



Two Users

- Ease of use; one-touch operation
- Extra-large LCD for easy reading of measurement
- Measurements consist of systolic and diastolic blood pressure and pulse rate
- Arrhythmia reminder
- Large memory space of 2 x 99 group of memories; for 2 users
- Performs average of last 3 readings
- 30s Auto shut off
- Bluetooth compatible to Jumper Health App

Specifications

JUMPER

Type	Arm Type
LCD Size (mm)	60 x 77
Cuff Length (cm)	22-36/ 22-42
User Memory	2 x 99 records
Automatic Shutdown	30s
Power Supply	AA x 4 USB 5V Adapter
Measurement Range	Pressure Value 0-295 ± 3mmHg Pulse Rate 40-199 ± 5% BPM



JUMPER Health App

Measure



FHR



Weight



Oximeter



Thermometer



Smart Wireless Thermometer



TENS Therapy



Blood pressure



Multi-Device Connectivity

Trend Graph

- Fast and Accurate data acquisition that reduces human error
- Easy Data Access
- Track Vital Signs without hassle

